

Introductory CyTOF®2 Operator Training Course


July 2013



Introductory CyTOF® 2 Operator Training Course

TABLE OF CONTENTS

- 1) Course Info
 - a) Course overview
 - b) Table of contents
 - c) Intro Training Feedback Form
- 2) Journey of Cell (PPT only)
- 3) Daily QC (PPT)
- 4) Operation Guide
- 5) Maintenance (PPT)
- 6) Manual Tuning (PPT)
- 7) Sample Acquisition (PPT)
- 8) Reagent Labeling
 - a) PPT
 - b) List of MAXPAR Kits
 - c) Kit Protocol printed
 - d) X8 and DN3 MAXPAR Kits Explanation
 - e) DVS Metal-conjugated Antibodies
 - f) Determining Metal Content for Antibody Conjugates
- 9) Special Topics
 - a) Sources of Signal Spillover (PPT)
 - b) Using Calibration Beads (PPT)
- 10) Protocols
 - a) CyTOF Experiment Planning Worksheet
 - b) CyTOF Surface-marker Staining Protocol Ms Splenocytes Draft_1.2
 - c) CyTOF Surface-marker Staining Protocol Hu PBMC Draft_1.4
 - d) Cisplatin Protocol_Viability MeOH Fix_Nolan Lab



DVS Sciences


CyTOF®2 Technical Overview

Utilizing the Power of the Atomic Spectrum

DVS Sciences

13 Lanthanides: 87 stable isotopes

24 elements: 67 stable isotopes



metal-tagged antibody

Lanthanides highlighted in green:


La	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
3	5	3	2	3	2	3	4	3	4	3	4	3	2

Samarium
6 isotopes

The detector has ninety three mass channels.
DVS sells 32 MAXPAR metal kits plus Rhodium and Iridium intercalator reagents.

Utilizing the Atomic Mass Spectrum


DVS Sciences



©Erin Simonds --

Atomic Mass Spectrum

DVS Sciences



- Abundant tags of similar intensity
- Discrete signals: minimal overlap
- Single metal controls not required
- Background cellular signal: zero

Mass Cytometry

DVS Sciences



CyTOP® 2 Mass Cytometer

- Inductively Coupled Plasma (ICP) Time-of-Flight (TOF) mass spectrometry discriminates metal-conjugated probes on a per-cell basis
- 103-193 Dalton atomic analytical mass range
- Single cell suspension input
- .fcs and .bt data output
- Optional Auto-sampler (96 well plates)




MAXPAR® Kits and Metal-conjugated antibodies

- Metal-conjugated antibodies (as of April 2013)
 - ~106 anti-human, 54 anti-mouse, 9 secondary antibodies
 - Reactive with cell surface, cytokine and signaling epitopes
- MAXPAR® Element labeling kits for antibody conjugations

Mass Cytometry Workflow

DVS Sciences



Bendall S, et al. (2012). Trends in Immunology, 33 (7), 323-332

Journey of the Cell: Outline

DVS Sciences

INPUT = Cells in liquid suspension stained with metal-conjugated probes


OUTPUT = Individual metal ions separated on the basis of mass

The CyTOF® 2 achieves this through the following steps:

1. Sample introduction and ionization
PURPOSE - to strip water from the cells followed by vaporization, atomization and ionization within the plasma
2. Ion filtering
PURPOSE - to filter out unwanted endogenous low mass ions and argon
3. Time of flight ion separation and detection of metal probes
PURPOSE - to separate the smallest from the highest mass ions; the time taken to reach the detector being proportional to mass

CyTOF® 2 Overview

DVS Sciences



Journey of the Cell: Sample Introduction and Ionization

DVS Sciences

INPUT
OUTPUT

The ICP-TOF-MS follows this through the following steps:

1. Sample introduction and ionization

PURPOSE - to strip water from the cells followed by vaporization, atomization and ionization within the plasma

Input: Single cell suspension

Output: Ions derived from metal-conjugated probes, endogenous cell components and argon


Sample introduction



PLASMA: to filter out unwanted endogenous ions, water molecules and argon

TOF: to separate the sample from the background ions based on the time taken to reach the detector being proportional to mass

Sample Introduction

DVS Sciences





Journey of the Cell: Plasma Generation

DVS Sciences


- Plasma is generated by inductively heating argon gas using an electrical coil.
- Argon gas flows through the torch, which itself is placed inside an induction coil supplied with a radio-frequency electric current.
- An electric spark is applied to the argon gas stream to induce free electrons, which interact with the radio-frequency magnetic field of the induction coil and are accelerated.
- The accelerated electrons collide with argon atoms, causing the release of an electron which is in turn accelerated by the magnetic field.
- This process continues, eventually converting the argon gas into high-temperature plasma consisting of argon atoms, free electrons and argon ions.



Inductively Coupled Plasma (ICP)

DVS Sciences

to vaporize, atomize and ionize the sample



Journey of the Cell: Ion Filtering

DVS Sciences

INPUT
OUTPUT

The DVS TORCH achieves this through the following steps:

INPUT: - Ionization from the gels followed by vaporization (desorption) and ionization in the plasma

Input

2. Ion filtering

PURPOSE - to filter out unwanted endogenous low mass ions and argon

Input: Ions derived from metal-conjugated probes, endogenous cell components and argon


Output: Highly enriched ions from the metal probes, randomly distributed in an ion cloud



Time of flight ion selection (TOF) is used for this purpose.

PURPOSE: To separate the ions derived from the tagged compounds, the low abundant ions from the detector during proportioning of the ions.

Ion Filtering


DVS Sciences





Journey of the Cell: Why do we Need a Vacuum? DVS Sciences

- The ions need to be in the gas-phase in order to freely pass through the system in response to the imposed electromagnetic fields
- A vacuum limits the number of random collisions and allows free movement of the ions
- A vacuum is sequentially applied as the ions travel through the cones, which removes ambient air from the system




Ion Filtering: Masses > 100 DVS Sciences

Quadrupole mass filter.



- Quadrupole Deflector:**
Eliminates non-ionized particles and photons
- Quadrupole Ion Filter:**
Eliminates ions with a mass below 100 Da
e.g. H⁺, C⁺, O⁺, N⁺, OH⁺, CO⁺, O₂⁺, Ar⁺, ArH⁺, ArO⁺ and Argon dimers



Journey of the Cell: TOF Ion Separation and Detection

DVS Sciences

3. Time of flight ion separation and detection of metal probes


PURPOSE - to separate the smallest from the highest mass ions; the time taken to reach the detector being proportional to mass



Input: Highly enriched ions from the metal probes, randomly distributed in an ion cloud



Output: Time-resolved ions from each individual metal probe

Time of Flight Ion Separation and Detection

DVS Sciences







Journey of the Cell: Summary

DVS Sciences

INPUT = Cells in liquid suspension stained with metal-conjugated probes

OUTPUT = Individual metal ions separated on the basis of mass

The CyTOF® 2 achieves this through the following steps:

1. Sample introduction and ionization

PURPOSE - to strip water from the cells followed by vaporization, atomization and ionization within the plasma

Output: Ions derived from metal-conjugated probes, endogenous cell components and argon

2. Ion filtering

PURPOSE - to filter out unwanted endogenous low mass ions and argon

Output: Highly enriched ions from the metal probes, randomly distributed in an ion cloud

3. Time of flight ion separation and detection of metal probes

PURPOSE - to separate the smallest from the highest mass ions; the time taken to reach the detector being proportional to mass

Output: Time-resolved ions from each individual metal probe