




CyTOF®2 Daily QC: Auto-Tuning

DVS Sciences

- Daily QC involves performance checking and tuning the CyTOF®2.
- Auto-Tuning automatically and simultaneously performs the following:
 - **Performance Check:** Testing instrument performance
 - **Tuning:** Calibration of CyTOF®2 for maximum performance
 - **Prior to QC:** Ensure CYTOF®2 is properly cleaned and assembled.
- Auto-Tuning is performed using **Tuning Solution**.

QC Product

DVS Sciences



CyTOF® Tuning Solution, E-Pure (250 ml), PN 201072 consists of six isotopes in liquid phase, each at the given concentration.

Isotopes Assessed in Auto-Tuning

DVS Sciences

Isotope:	Source:
^{133}Cs	Tuning Solution/Environmental contaminant
^{139}La	Tuning solution
^{155}Gd	Oxide formed from ^{139}La
^{159}Tb	Tuning solution
^{169}Tm	Tuning solution
^{193}Ir	Tuning solution

Auto-Tuning Results

DVS Sciences

Tuning: GC Results

Profiles: General Parameters | Masses Parameters | Dual Pulse Parameters | Detector V Parameters | Gases and Carrier Parameters | XE Parameters | GC Parameters | Control | Print All

Masses:

Mass - Time of Flight (M): 842.704

Mass - Time of Flight (T): 151.907

Resolution: 500

Instrument:

Optimal Detector V: 0.090

Optimal Nebulizer Gas: 0.25

Optimal Makeup Gas: 0.51

Optimal Carrier: 7.0

Optimal X: 0

Optimal Y: 0

Default: M11.M2: 0.0208

Analyte Name	Symbol	Mass	Dual Slope	Dual Intercept
Ca-40	Ca	132.905	0.02085	0.00000
La-138	La	138.905	0.02091	0.00000
Tl-153	Tl	152.925	0.02090	0.00000
Tl-163	Tl	162.934	0.02210	0.00000
Bi-190	Bi	190.962	0.02217	0.00000

Analyte	Symbol	Mass	Mean Intensity	Mean Duty	Mean Pulses	RSD Intensity %	RSD Pulses %	RSD Duty %
Ca-40	Ca	132.905	10,715,362	579,276	361,585	0.33	0.19	0.33
La-138	La	138.905	15,000,290	474,237	34,176	0.16	0.52	0.15
Tl-153	Tl	152.925	1,108,808	34,385	30,562	2.82	2.13	2.82
Tl-163	Tl	162.934	21,200,912	1,145,952	505,605	0.65	2.32	0.60
Bi-190	Bi	190.962	12,581,140	1,045,020	590,448	0.65	0.40	0.65
Bi-190	Bi	190.962	8,608	567	561	0.46	2.36	2.36
Bi-190	Bi	190.962	1,050,135	220,220	200,52	0.74	0.45	0.70

DVS Sciences

CyTOF[®]2 Auto-Tuning

CyTOF® 2 Daily QC: Auto-Tuning

DVS Sciences


- The purpose of QC is to check and optimize the following parameters:
 - Mass Calibration (Automatically run)
 - Mass Resolution (Automatically run)
 - Detector Voltage (selectable)
 - Dual Pulse Calibration (Automatically run)
 - X-Y Alignment (selectable)
 - Gases and Current Optimization (selectable)
 - QC Report (selectable)

Metal Signal Measurement


DVS Sciences

- Two different measurements are recorded for each metal in a 13 μ sec push:
 - Pulse count
 - Most accurate when ion counts are low enough that discrete ion pulses can be counted.
 - Detects the exact number of ion pulses per push.
 - Intensity
 - Used when pulse count is saturated, i.e. when ion count is high.
 - Detects the total intensity in the channel by measuring the area under the curve.

Pulse count:



Intensity:



Check Mass Resolution

DVS Sciences


- Mass Resolution

PURPOSE: To ensure that ions are tightly segregated by mass. Intensity curves should NOT overlap.

Time of Flight Separation of Ions

DVS Sciences


All ions of a particular mass reach the detector over a time range of 21 to 26 ns. This is known as the **integration time**. It is metal-specific.



Integration Time: metal specific

DVS Sciences

Symbol	Mass	Integration Time
Xe	131	21
Ce	133	21
Nd	142	22
Sm	147	22
Eu	151	23
Eu	153	23
Sm	154	23
Tb	159	23
Ho	165	24
Er	167	24
Yb	171	24
Yb	172	24




ToF values are in ns.

Check Mass Resolution

DVS Sciences

$$\text{Mass Resolution} = \frac{\text{ToF at } 159\text{Tb Max Intensity}}{2 \times \text{FWHM}}$$



Mass Resolution
above 400 = PASS

Check Mass Calibration


DVS Sciences

- Mass Calibration

PURPOSE: To ensure that ions are being detected in the correct detection channels on the TOF scale.

Check Mass Calibration

DVS Sciences




- The TOF for a particular ion is dependent on its mass to charge ratio (M/z). Because all ions have the same charge, the TOF for each ion is directly proportional to mass.
- The TOF value for each mass is **calibrated at install** to ensure optimal measurement in the appropriate detector channel.
- Auto-Tuning checks Time of Flight for the ^{133}Cs and ^{193}Ir peaks and then calculates values for additional isotopes.

PURPOSE: To ensure that ions are being detected in the correct detection channel on the TOF scale.

Check Mass Calibration

DVS Sciences



Element	Mass (Da)	TOF (peak) (nsec)
Cs	133	9,533
La	139	9,786
Tb	159	10,591
Tm	169	10,975
Ir	191	11,781
Ir	193	11,852

Detector Voltage Optimization

DVS Sciences

- Detector Voltage

PURPOSE: To ensure that the system is optimized for detector sensitivity and longevity.

Detector Voltage Optimization

DVS Sciences


- The detector must be optimized for the correct ratio between pulse count and intensity; this is important for instrument dual pulse calibration, absolute quantification and dynamic range.
- The detector must also be powered with the correct number of volts: this will increase as the detector ages and sensitivity decreases.

Detector Voltage Optimization

DVS Sciences

Dual Slope Stability Method

- Intensity is determined for one, two and three ion counts at Detector V1 and Detector V2.
- The optimal voltage is then interpolated to give the target dual slope of 0.03.



Dual Count Calibration

DVS Sciences

- Dual Pulse Calibration

Purpose: To determine a Dual Pulse Coefficient (Dual Slope) that correlates Pulse Count and Intensity.

What are Dual Counts?

DVS Sciences

- The TOF detector sees a wide range of signal (low to high analyte density)
- **Pulse Count** (for low concentrations) and **Intensity** (for high concentrations) values are collected for every channel.
- CyTOF plots the entire data range on a single Dual Signal scale, the *units of which are actual counts of particles that hit the detector*. This requires:
 - Determination of a Dual Count Coefficient which relates Pulse Count and Intensity values.
 - Application of the following equation to determine actual counts:

$$\text{Counts} = \text{Intensity} \times \text{Dual Count Coefficient}$$

- The instrument-derived (Di) Dual Count Coefficient is determined using the Dual Pulse Calibration during Auto-Tuning of the instrument.

Dual Pulse Calibration

DVS Sciences




Figure 3: Derivation of the Dual Count Coefficient using Instrument Dual Count Calibration

During the calibration, tandem solution dual count data are collected across a wide range of quadrupole deflection voltages, from least to most optimal (left to right). At low optimal deflection, analyte count rates tend to be low and stable data (Pulse Count - Intensity $R^2 < 0.7$). At optimal deflection, analyte count rates causes pulse count saturation. In between is a zone of dual validity in which the pulse count and intensity are linearly related. The slope of the line in this region is the Dual Count Coefficient.

- Dual Count Coefficients (DCC) for two analytes at opposite ends of the mass range (^{133}Cs and ^{169}Tm) are determined.
- An equation is determined relating mass and DCC.
- The DCC for all other masses are extrapolated from this equation.

X-Y Alignment

DVS Sciences


- X-Y Alignment

PURPOSE: To fine tune the position of the interface to ensure that ion clouds are being optimally pulled into the cones.

X-Y Alignment

DVS Sciences

- Auto-Tuning runs an algorithm which uses spiral cycles that begin from a pre-determined starting point.
- During each spiral cycle, the software determines the point where the highest ^{159}Tb signal is achieved.
- The next spiral cycle begins from the point where the high ^{159}Tb signal was found.



Gases and Current Optimization

DVS Sciences

- Gases and Current Optimization

PURPOSE: To ensure that plasma temperature is high enough and that conditions are optimal for ion clouds to travel into the cones.

Gases Optimization

DVS Sciences


- The plasma temperature is 5500° K.
- This must be maintained for optimum vaporization, atomization and ionization.
- The formation of metal oxides ($M + 16$) increases when the plasma temperature drops.

PURPOSE: To ensure that the plasma temperature is high enough.

Gases Optimization

DVS Sciences

- A percentage of Lanthanum oxide (measured in the 155 channel) above 3% indicates the plasma temperature is too low.



PURPOSE: To ensure that the plasma temperature is high enough.

Current Optimization

DVS Sciences

- Current is checked at optimized Make-up Gas and Nebulizer Gas values.
- The optimum current is where the maximum ^{159}Tb signal is achieved and oxides are still below 3%.

PURPOSE: To ensure that the current on cones is optimal to drive ion clouds through the interface .

QC Report

DVS Sciences

- QC Report – ^{159}Tb Dual Mean and RSD %*

PURPOSE: To ensure that CyTOF® 2 is optimally detecting ions .

* RSD = Relative Standard Deviation (equivalent to CV)